

Product Datasheet

Mouse PCSK9 ELISA Kit (orb441449)

Description	Mouse Proprotein Convertase Subtilisin/Kexin Type 9 (PCSK9) Sandwich ELISA Kit
Reactivity	Mouse
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q80W65
Sensitivity	6.4pg/mL
Entrez	100102
Expiration Date	Please enquire.