

Product Datasheet

Human PDCD1LG1 ELISA Kit (orb441459)

Description	Human Programmed Cell Death Protein 1 Ligand 1 (PDCD1LG1) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q9NZQ7
Sensitivity	0.066ng/mL
Entrez	29126
Expiration Date	Please enquire.