

Product Datasheet

Human PDCD1LG2 ELISA Kit (orb441461)

Description	Human Programmed Cell Death Protein 1 Ligand 2 (PDCD1LG2) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9BQ51
Sensitivity	0.063ng/mL
Entrez	80380
Expiration Date	Please enquire.