

Product Datasheet

Rat PDCN ELISA Kit (orb441466)

Description	Rat Podocin (PDCN) Sandwich ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	urine, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	urine, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q8K4G9
Sensitivity	0.053ng/mL
Entrez	170672
Expiration Date	Please enquire.