

Product Datasheet

Human PDHa ELISA Kit (orb441473)

Description	Human Pyruvate Dehydrogenase Alpha (PDHa) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P08559
Sensitivity	0.078ng/mL
Entrez	5160
Expiration Date	Please enquire.