

Product Datasheet

Human PDIA4 ELISA Kit (orb441481)

Description	Human Protein Disulfide Isomerase A4 (PDIA4) Sandwich ELISA Kit
Reactivity	Human
Range	1.56-100ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P13667
Sensitivity	0.66ng/mL
Entrez	9601
Expiration Date	Please enquire.