

Product Datasheet

Mouse PDPK1 ELISA Kit (orb441489)

Description	Mouse Phosphoinositide Dependent Protein Kinase 1 (PDPK1) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9Z2A0
Sensitivity	0.055ng/mL
Entrez	18607
Expiration Date	Please enquire.