

Product Datasheet

Human PECAM1 ELISA Kit (orb441496)

Description	Human Platelet/Endothelial Cell Adhesion Molecule (PECAM1) Sandwich ELISA Kit
Reactivity	Human
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P16284
Sensitivity	13.6pg/mL
Entrez	5175
Expiration Date	Please enquire.