

Product Datasheet

Mouse PECAM1 ELISA Kit (orb441497)

Description	Mouse Platelet/Endothelial Cell Adhesion Molecule (PECAM1) Sandwich ELISA Kit
Reactivity	Mouse
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	plasma.
Sample Types	plasma.
Uniprot ID	Q08481
Sensitivity	6.5pg/mL
Entrez	18613
Expiration Date	Please enquire.