

Product Datasheet

Rat PER2 ELISA Kit (orb441507)

Description	Rat Period Circadian Protein 2 (PER2) Sandwich ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q9Z301
Sensitivity	0.101ng/mL
Entrez	63840
Expiration Date	Please enquire.