

Product Datasheet

Human PF4V1 ELISA Kit (orb441516)

Description	Human Platelet Factor 4 Variant 1 (PF4V1) Sandwich ELISA Kit
Reactivity	Human
Range	1.56-100ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P10720
Sensitivity	0.8ng/mL
Entrez	5197
Expiration Date	Please enquire.