

Product Datasheet

Human PFKM ELISA Kit (orb441519)

Description	Human Phosphofructokinase, Muscle (PFKM) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P08237
Sensitivity	0.114ng/mL
Entrez	5213
Expiration Date	Please enquire.