

Product Datasheet

Chicken PGC ELISA Kit (orb441527)

Description	Chicken Pepsinogen C (PGC) Sandwich ELISA Kit
Reactivity	Gallus
Range	0.781-50ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q9PWK1
Sensitivity	0.198ng/mL
Entrez	395690
Expiration Date	Please enquire.