

Product Datasheet

Human PGD2S ELISA Kit (orb441531)

Description	Human Prostaglandin D2 Synthase, Brain (PGD2S) Sandwich ELISA Kit
Reactivity	Human
Range	62.5-4000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, cerebrospinal fluid and other biological fluids.
Sample Types	serum, plasma, cerebrospinal fluid and other biological fluids.
Uniprot ID	A0A024R8G3
Sensitivity	26.5pg/mL
Entrez	5730
Expiration Date	Please enquire.