

Product Datasheet

Human PGLYRP1 ELISA Kit (orb441539)

Description	Human Peptidoglycan Recognition Protein 1 (PGLYRP1) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, cell lysates and other biological fluids.
Sample Types	serum, plasma, cell lysates and other biological fluids.
Uniprot ID	075594
Sensitivity	0.115ng/mL
Entrez	8993
Expiration Date	Please enquire.