

Product Datasheet

Rat PHB ELISA Kit (orb441548)

Description	Rat Prohibitin (PHB) Sandwich ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P67779
Sensitivity	0.121ng/mL
Entrez	25344
Expiration Date	Please enquire.