

Product Datasheet

Human PIIINP ELISA Kit (orb441565)

Description	Human Procollagen III N-Terminal Propeptide (PIIINP) Sandwich ELISA Kit
Reactivity	Human
Range	15.6-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum and other biological fluids.
Sample Types	serum and other biological fluids.
Sensitivity	5.60pg/mL
Expiration Date	Please enquire.