

Product Datasheet

Human PKCz ELISA Kit (orb441603)

Description	Human Protein Kinase C Zeta (PKCz) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q05513
Sensitivity	0.112ng/mL
Entrez	5590
Expiration Date	Please enquire.