

Product Datasheet

Mouse PLA2G2A ELISA Kit (orb441623)

Description	Mouse Phospholipase A2, Group IIA (PLA2G2A) Sandwich ELISA Kit
Reactivity	Mouse
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P31482
Sensitivity	5.7pg/mL
Entrez	18780
Expiration Date	Please enquire.