

Product Datasheet

Human PLA2G4D ELISA Kit (orb441628)

Description	Human Phospholipase A2, Group IVD (PLA2G4D) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q86XP0
Sensitivity	0.064ng/mL
Entrez	280748
Expiration Date	Please enquire.