

Product Datasheet

Human PLAA ELISA Kit (orb441631)

Description	Human Phospholipase A2 Activating Protein (PLAA) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q9Y263
Sensitivity	0.114ng/mL
Entrez	9373
Expiration Date	Please enquire.