

Product Datasheet

Human PLB ELISA Kit (orb441634)

Description	Human Phospholipase B (PLB) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q6P1J6
Sensitivity	0.113ng/mL
Entrez	151056
Expiration Date	Please enquire.