

Product Datasheet

Human PLD2 ELISA Kit (orb441638)

Description	Human Phospholipase D2 (PLD2) Sandwich ELISA Kit
Reactivity	Human
Range	1.56-100ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	O14939
Sensitivity	0.58ng/mL
Entrez	5338
Expiration Date	Please enquire.