

Product Datasheet

Human CSH1 ELISA Kit (orb441643)

Description	Human Chorionic Somatomammotropin Hormone 1 (CSH1) Sandwich ELISA Kit
Reactivity	Human
Range	3.125-200ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	PODML2
Sensitivity	1.12ng/mL
Entrez	1442
Expiration Date	Please enquire.