

Product Datasheet

CD56 antibody (orb44165)

Description

Mouse Monoclonal to CD56.

Reactivity

Human, Primate

Conjugation

Unconjugated

Tested Applications

FC, ICC, IHC-P, IP

Immunogen

KG-1 human acute myelogenous leukemia cell line

Target

CD56

Preservatives

Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide

Concentration

1 mg/ml

Storage

Store at 2-8°C. Do not freeze.

Note

For research use only

Application notes

Flow cytometry: Recommended dilution: 3 µg/ml.
Immunohistochemistry (paraffin sections): Recommended dilution: 1:50; pretreatment: heat retrieval (microwave oven) in citrate buffer, standard ABC visualization; positive tissue: membrane-specific reaction in adjacent peripheral nerve; staining of paraprostatic autonomic ganglion.

Isotype

Mouse IgG2a

Clonality

Monoclonal

Purity

Purified by protein-A affinity chromatography.

Uniprot ID

P13591

Entrez


4684

Dilution Range

Flow cytometry: Recommended dilution: 3 µg/ml.
Immunohistochemistry (paraffin sections): Recommended dilution: 1:50; pretreatment: heat retrieval (microwave oven) in citrate buffer, standard ABC visualization; positive tissue: membrane-specific reaction in adjacent peripheral nerve; staining of paraprostatic autonomic ganglion.

Expiration Date

12 months from date of receipt.



Flow
cytometry
analysis
(surface
stainin...