

Product Datasheet

Mouse Lipase, Pancreatic ELISA Kit (orb441650)

Description	Mouse Lipase, Pancreatic (PL) Sandwich ELISA Kit
Reactivity	Mouse
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q6P8U6
Sensitivity	35pg/mL
Entrez	69060
Expiration Date	Please enquire.