

Product Datasheet

Human PLOD3 ELISA Kit (orb441653)

Description	Human Procollagen Lysine-2-Oxoglutarate-5-Dioxygenase 3 (PLOD3) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	O60568
Sensitivity	0.106ng/mL
Entrez	8985
Expiration Date	Please enquire.