

Product Datasheet

Rat Lipase, Pancreatic ELISA Kit (orb441657)

Description	Rat Lipase, Pancreatic (PL) Sandwich ELISA Kit
Reactivity	Rat
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P27657
Sensitivity	37pg/mL
Entrez	25702
Expiration Date	Please enquire.