

Product Datasheet

Human PPARα ELISA Kit (orb441692)

Description	Human Peroxisome Proliferator Activated Receptor Alpha (PPARα) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q07869
Sensitivity	0.052ng/mL
Entrez	5465
Expiration Date	Please enquire.