

Product Datasheet

Human PPAR γ C1a ELISA Kit (orb441698)

Description	Human Peroxisome Proliferator Activated Receptor Gamma Coactivator 1 Alpha (PPAR γ C1a) Sandwich ELISA...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9UBK2
Sensitivity	0.061ng/mL
Entrez	10891
Expiration Date	Please enquire.