

Product Datasheet

Porcine pPLA2 ELISA Kit (orb441707)

Description	Porcine Phospholipase A2, Pancreas (pPLA2) Sandwich ELISA Kit
Reactivity	Porcine
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P00592
Sensitivity	0.059ng/mL
Entrez	445525
Expiration Date	Please enquire.