

Product Datasheet

Mouse PPL ELISA Kit (orb441709)

Description	Mouse Periplakin (PPL) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9R269
Sensitivity	0.115ng/mL
Entrez	19041
Expiration Date	Please enquire.