

Product Datasheet

Human PPM1A ELISA Kit (orb441710)

Description	Human Protein Phosphatase, Mg ²⁺ /Mn ²⁺ Dependent 1A (PPM1A) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P35813
Sensitivity	0.118ng/mL
Entrez	5494
Expiration Date	Please enquire.