

Product Datasheet

Rat PPP2R4 ELISA Kit (orb441715)

Description	Rat Protein Phosphatase 2A Activator, Regulatory Subunit 4 (PPP2R4) Sandwich ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	B2RYQ2
Sensitivity	0.057ng/mL
Entrez	362102
Expiration Date	Please enquire.