

Product Datasheet

Human PRDX5 ELISA Kit (orb441733)

Description	Human Peroxiredoxin 5 (PRDX5) Sandwich ELISA Kit
Reactivity	Human
Range	0.781-50ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cerebrospinal fluid and other biological fluids.
Sample Types	tissue homogenates, cerebrospinal fluid and other biological fluids.
Uniprot ID	P30044
Sensitivity	0.34ng/mL
Entrez	25824
Expiration Date	Please enquire.