

Product Datasheet

Human PRKAb1 ELISA Kit (orb441743)

Description	Human Protein Kinase, AMP Activated Beta 1 (PRKAb1) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q9Y478
Sensitivity	0.124ng/mL
Entrez	5564
Expiration Date	Please enquire.