

Product Datasheet

Rat PRKAb1 ELISA Kit (orb441744)

Description	Rat Protein Kinase, AMP Activated Beta 1 (PRKAb1) Sandwich ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P80386
Sensitivity	0.060ng/mL
Entrez	83803
Expiration Date	Please enquire.