

Product Datasheet

Human PRKG1 ELISA Kit (orb441746)

Description	Human Protein Kinase, cGMP Dependent Type I (PRKG1) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q13976
Sensitivity	0.063ng/mL
Entrez	5592
Expiration Date	Please enquire.