

Product Datasheet

Bovine PRL ELISA Kit (orb441748)

Description	Bovine Prolactin (PRL) Sandwich ELISA Kit
Reactivity	Bovine
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P01239
Sensitivity	0.066ng/mL
Entrez	280901
Expiration Date	Please enquire.