

Product Datasheet

Rat PROC ELISA Kit (orb441760)

Description	Rat Protein C (PROC) Sandwich ELISA Kit
Reactivity	Rat
Range	1.56-100ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	plasma.
Sample Types	plasma.
Uniprot ID	P31394
Sensitivity	0.42ng/mL
Entrez	25268
Expiration Date	Please enquire.