

Product Datasheet

Human PSA ELISA Kit (orb441782)

Description	Human Puromycin Sensitive Aminopeptidase (PSA) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P55786
Sensitivity	0.118ng/mL
Entrez	9520
Expiration Date	Please enquire.