

Product Datasheet

Rat PSGL1 ELISA Kit (orb441790)

Description	Rat P-Selectin Glycoprotein Ligand 1 (PSGL1) Sandwich ELISA Kit
Reactivity	Rat
Range	0.625-40ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q8K5B0
Sensitivity	0.227ng/mL
Entrez	363930
Expiration Date	Please enquire.