

Product Datasheet

Human PSMD10 ELISA Kit (orb441795)

Description	Human Proteasome 26S Subunit, Non ATPase 10 (PSMD10) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	075832
Sensitivity	0.059ng/mL
Entrez	5716
Expiration Date	Please enquire.