

Product Datasheet

Canine PTH ELISA Kit (orb441819)

Description	Canine Parathyroid Hormone (PTH) competitive inhibition ELISA Kit
Reactivity	Canine
Range	12.5-800pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P52212
Sensitivity	3.52pg/mL
Entrez	403986
Expiration Date	Please enquire.