

Product Datasheet

Human PTPLA ELISA Kit (orb441834)

Description	Human Protein Tyrosine Phosphatase Like Protein A (PTPLA) Sandwich ELISA Kit
Reactivity	Human
Range	0.781-50ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	B0YJ81
Sensitivity	0.25ng/mL
Entrez	9200
Expiration Date	Please enquire.