

Product Datasheet

Human PTPN6 ELISA Kit (orb441837)

Description	Human Protein Tyrosine Phosphatase, Non Receptor Type 6 (PTPN6) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P29350
Sensitivity	0.105ng/mL
Entrez	5777
Expiration Date	Please enquire.