

Product Datasheet

Human PTX3 ELISA Kit (orb441842)

Description	Human Pentraxin 3, Long (PTX3) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, saliva and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, saliva and other biological fluids.
Uniprot ID	P26022
Sensitivity	0.114ng/mL
Entrez	5806
Expiration Date	Please enquire.