

Product Datasheet

Rat PUMA ELISA Kit (orb441845)

Description	Rat p53 Upregulated Modulator Of Apoptosis (PUMA) Sandwich ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q80ZG6
Sensitivity	0.047ng/mL
Entrez	317673
Expiration Date	Please enquire.