

Product Datasheet

Human PVRL4 ELISA Kit (orb441850)

Description	Human Poliovirus Receptor Related Protein 4 (PVRL4) Sandwich ELISA Kit
Reactivity	Human
Range	125-8000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q96NY8
Sensitivity	52pg/mL
Entrez	81607
Expiration Date	Please enquire.