

Product Datasheet

Human PXR ELISA Kit (orb441852)

Description	Human Pregnane X Receptor (PXR) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	075469
Sensitivity	0.126ng/mL
Entrez	8856
Expiration Date	Please enquire.